



# Thermal Decomposition Graphs – TGA & DSC (16 Marks)

---

## ◆ 1. Introduction

- **Thermal Analysis** is used to study how materials behave when heated.
  - Two key methods:
    - **TGA (Thermogravimetric Analysis)**: Measures **weight change** with temperature.
    - **DSC (Differential Scanning Calorimetry)**: Measures **heat flow** during heating/cooling.
- 

## ◆ 2. TGA – Thermogravimetric Analysis

### Principle:

- Measures **mass loss** of a material as temperature increases.
- Used to detect **decomposition, oxidation, and evaporation**.

### Graph:

- **X-axis**: Temperature (°C or K)
- **Y-axis**: % Weight (mass)
- The curve shows **mass decreasing** with temperature.

### Stages:

1. **Initial flat region**: No change – stable.
2. **Drop region**: Sudden weight loss – decomposition or evaporation.
3. **Plateau**: No further reaction – residue remains.

### Applications:

- Checking **thermal stability**.
  - Studying **decomposition temperature**.
  - Measuring **moisture, volatile content, and residue**.
- 

## ◆ 3. DSC – Differential Scanning Calorimetry

### Principle:

- Measures **heat flow difference** between sample and reference as they are heated/cooled.
- Identifies **endothermic** (heat absorbed) or **exothermic** (heat released) events.

### Graph:

- **X-axis**: Temperature (°C)

- **Y-axis:** Heat Flow (mW or mJ/sec)
- Peaks indicate thermal transitions.

#### **Types of Peaks:**

- **Endothermic peak (downward):** Melting, dehydration, glass transition.
- **Exothermic peak (upward):** Crystallization, curing, combustion.

#### **Applications:**

- **Melting point, crystallization temperature.**
- Studying **glass transition**.
- **Polymer and drug** analysis.

### 4. Diagram Description

#### **TGA Graph:**

- Straight line → Sudden drop → Flat line.
- Label: "Mass (%) vs Temperature (°C)"
- Drop = decomposition or weight loss.

#### **DSC Graph:**

- Baseline → Peak (up/down) → Return to baseline.
- Label: "Heat Flow vs Temperature"
- Downward peak = melting (endothermic)
- Upward peak = crystallization (exothermic)

### 5. Advantages

TGA	DSC
Simple setup	Precise heat flow measurement
Measures mass change	Identifies phase transitions
Useful for residue and moisture	Determines melting & glass points

### 6. Disadvantages

TGA	DSC
Cannot identify chemical changes	No weight change is detected
Needs inert gas for accuracy	Sensitive to sample mass & heating rate

## 8. Summary

- **TGA** shows **weight loss** vs temperature → used for decomposition and stability study.
  - **DSC** shows **heat flow** vs temperature → used for melting, glass transition, and curing.
  - Both are essential in **thermal analysis** of materials.
-